

## East SP & Zeus SP Flexible Pendants



*The SP Flexible Pendant range provides a cost-effective solution to trailing hoses around the operating table*

- ◆ Complies with BS EN 737-1, BS 5682:1984, HTM 2022 and C11
- ◆ Up to six outlets available
- ◆ Colour coded hoses and terminal units
- ◆ Integral first fix check valves
- ◆ Non-interchangeable hose connections
- ◆ User-definable hose length
- ◆ East SP - all metal construction
- ◆ Zeus SP - metal first fix, polymer second fix
- ◆ Three year warranty

Partnership for Life

## Product Summary

The Penlon SP Flexible Pendant shall comprise colour coded flexible hoses, with medical service terminal units, suspended from gas specific first fix assemblies. Terminal units and first fix assemblies will be protected by white decorative covers.

The terminal units shall comply with BS EN 737 and accept probes to BS 5682:1984 or BS 6834:1987(1992). The terminal units will not incorporate anti-swivel pins.

## Product Features

### Design features

- ◆ Meets the requirements of BS EN 737-1, BS 5682:1984, HTM 2022 and C11
- ◆ CE marked
- ◆ Colour coded terminal units
- ◆ Anti-static, colour coded, reinforced hoses
- ◆ Hose length to customer requirements
- ◆ Hose and terminal unit assemblies can be removed without isolating pipeline

### Easy to install

- ◆ Available as a complete kit, or as separate first and second fix components so that supply can be matched to site requirements
- ◆ First fix connectors incorporate copper stub pipes for ease of installation using inert gas jointing procedures
- ◆ Second fix comprises of hose and terminal unit assemblies

Colour coded hose assemblies will be manufactured with non-interchangeable male and female terminations permanently attached. Hose lengths will be suited to specific site requirements.

First fix connectors shall be manufactured from machined brass, with copper stub pipes for easy installation, and will incorporate check valves to allow the removal of hoses without the need to disrupt the gas supply.

### Gases and Services

- ◆ Oxygen
- ◆ Nitrous Oxide
- ◆ 50% Oxygen / 50% Nitrous Oxide
- ◆ Medical Air (400 kPa)
- ◆ Surgical Air (700 kPa)
- ◆ Medical Vacuum
- ◆ AGSS
- ◆ Nitrogen Contact Sales Office for further information
- ◆ Carbon Dioxide Contact Sales Office for further information

### Service and Maintenance

- ◆ Three year warranty from date of manufacture, excluding components designated as service items to be replaced by the customer as part of the published service schedule
- ◆ Throughout the UK our network of service engineers provides full Service and Maintenance support of the highest quality
- ◆ Outside the UK our fully trained distributors provide a similar local service

### Component List (complete kit)

- ◆ Gas specific first fix ceiling NIST assemblies
- ◆ Colour coded hoses and terminal units
- ◆ Installation kit

## Ordering Information

Contact our Sales Office, or refer to our latest price list for part number and price information

### Medical Gas Pipeline Sales

Tel: +44 (0) 1235 547038

Fax: +44 (0) 1235 547055

E-Mail: [medgas@penlon.co.uk](mailto:medgas@penlon.co.uk)

### Website

<http://www.penlon.com/pipeline.html>



Penlon Limited, Abingdon Science Park,  
Barton Lane, Abingdon OX14 3PH, UK

Penlon Limited reserves the right to change technical specification without prior notice

© Penlon Limited 2004

