

## CAVS and R-8000 Resuscitation Carts



*The CAVS and R-8000 provide a choice of well designed, versatile life support carts that can be equipped to the specific requirements of the user*

- ◆ CAVS Cart
  - Compact and functional, minimum footprint, maximum storage capacity
  - Multi-drawer unit
  - Lightweight but rigid polypropylene construction
- ◆ R-8000 Resuscitation Cart
  - Top rail plus two side rails provides multi-option mounting capacity
  - Storage for two cylinders
  - Deep, high capacity drawers with secure locks

Partnership for Life

## R-8000 Resuscitation Cart

A well-engineered, high quality resuscitation cart

You can load the R-8000 with all the equipment specific to your requirements, and hold a stock of drugs in its secure, lockable draw units. You then have a mobile unit ready for instant response to emergency situations. The R-8000 is easy to manoeuvre; noiseless and stable, and built to last.

- ◆ Steel construction, coated with epoxy resin, mounted on non-skid castors, 125 mm diameter. Two braked, plus one antistatic
- ◆ Two drawers (310 x 400 x 470 mm) with dividers and keys, additional two drawers available as an option
- ◆ Tray (450 x 400 mm), supports up to 30 kg
- ◆ Two side rails, plus one top rail for gas therapy and suction equipment: (e.g. oxygen regulator, flowmeter, humidifier, suction controller and collection jar)
- ◆ Secure support for two cylinders
- ◆ Optional cardiac massage board
- ◆ Drip stand, which can be positioned at either side of the cart



R-8000 with optional four drawer unit and cardiac massage board, other devices must be ordered separately



**Castors**  
Four non-skid, 125 mm diameter, two braked and one antistatic



**Cardiac massage board**  
The optional cardiac massage board can be stored on the side of the cart



**High level rail**  
Standard configuration rail for oxygen therapy and suction devices (available separately)



**Top tray**  
Heavy duty tray supports up to 30 kg



**Drawer units**  
Two lockable drawer units (310 x 400 x 470 mm) with dividers and keys

# CAVS

## Cart for Advanced Vital Support

A modern, multipurpose cart, for theatres and wards

In theatre, the seven large drawer units make CAVS an ideal mobile storage unit for sterile instruments for intubation and surgery. For general emergency situations CAVS can be loaded with a defibrillator, oxygen cylinder, and drip system.

- ◆ Small footprint takes up minimal space in your work area
- ◆ Smooth, easy to manoeuvre operation
- ◆ Sturdy upper tray for defibrillator equipment
- ◆ MRI safe
- ◆ Secure mount for oxygen cylinder, and parallel fittings for suction systems
- ◆ Drawer unit with seven drawers, adjustable dividers and labels
- ◆ Left- or right-hand side drip stand mounting, plus low-level storage point



**Cylinder storage**  
Secure location for one size D/E oxygen cylinder



**Drawers**  
Seven draw units, with adjustable dividers and labels



**Top tray**  
Heavy duty tray and left- or right-hand side mounting for drip stand



**Side tray**  
Folding side tray gives added work space if required



**Castors**  
Four non-skid, 100 mm diameter castors, two braked, and one antistatic

## Technical Specification

	R-8000	CAVS
Size (H x W x D)	1400 x 470 x 680 mm	930 x 430 x 570 mm
Weight	56 kg	10 kg
Construction	Steel, coated with epoxy resin	Polypropylene
Castors	Four, non-skid 125 mm diameter castors, two braked and one anti-static	Four, non-skid 100 mm diameter castors, two braked and one anti-static
Drawers	Two lockable drawers with dividers and keys. Additional two-drawer unit available as an option	Seven, with adjustable dividers and labels
Medical Rail	Two side and one high level rail sections	N/A
Cylinder Storage	2 x cylinders (not supplied)	1 x size D or E cylinder (not supplied)
Cardiac Massage Board	Optional	No
Drip Stand	Left or right-hand side mounting, supplied	Left or right-hand side mounting, supplied
MRI Safe	No	Yes

## Ordering Information

### R-8000

Part No.	Description
4230010	Resuscitation Cart (includes two drawers and drip stand)
3230014	Additional two-drawer unit
0200020	Cardiac massage board
1761-000	Rail-mounted mains electrical four-way adaptor
1851-000	Rail-mounted rubbish bag holder

### CAVS

Part No.	Description
5230505	CAVS Cart (includes seven drawers, labels, and drip stand)
4230505	Mounting plate for defibrillator
4230507	Support straps for use with mounting plate
4230504	Pull-out tray
4230503	Folding side tray
4230508	Support bracket for compact defibrillator / monitor

For suction and oxygen therapy products, and additional rail-mounted accessories, please contact the Sales Office at Penlon Limited

Information contained in this leaflet is correct at the date of publication. The policy of Penlon Limited is one of continued improvement to its products. Because of this policy Penlon Limited reserves the right to make any changes, which may affect the information in this leaflet without giving prior notice.

Distributed by: Penlon Limited • Manufactured by: Hersill, S.L.

[www.penlon.com](http://www.penlon.com)

Penlon Limited Abingdon Science Park Barton Lane Abingdon OX14 3PH UK

#### General

Tel: +44 (0) 1235 547000  
Fax: +44 (0) 1235 547041  
<http://www.penlon.com>

#### International Sales

Tel: +44 (0) 1235 547001  
Fax: +44 (0) 1235 547021  
E-Mail: [export@penlon.co.uk](mailto:export@penlon.co.uk)

#### UK Sales

Tel: 01235 547036  
Fax: 01235 547023  
E-Mail: [uksales@penlon.co.uk](mailto:uksales@penlon.co.uk)

#### Service

Tel: 01235 547060  
Fax: 01235 547061  
E-Mail: [service@penlon.co.uk](mailto:service@penlon.co.uk)

#### Medical Gas

Tel: 01235 547038  
Fax: 01235 547055  
E-Mail: [medgas@penlon.co.uk](mailto:medgas@penlon.co.uk)